

Sentida SC Bed



**Customisable & efficient
nursing care bed for
demanding environments.**

Sentida SC Bed



New concepts for living in old age are in demand. These include the desire to be provided with more homeliness, individuality, safety, comfort and flexibility in the equipment of corresponding living arrangements.

With the new nursing bed platform sentida sc and customied concepts, we can respond to the different requirements even more flexibly, efficiently and demand-oriented. We want to support the maintenance and promotion of individual mobility as an overarching and essential goal in care in a solution-oriented manner.

Product Benefits:

- Fully customisable with various models, decor and designs
- SafeFree® side rails – safe mobilization regardless of body size.
- MobiStick 2 – manual mobilization assistant
- SafeControl® handset – clear, ergonomic and intuitive
- Double retraction of the lying surface – safe and comfortable mobilisation in bed
- Automatic positions – quick, easy and ergonomic thanks to the one-button philosophy.
- SafeFree® in orientation position – fear reduction without unnecessary barriers.
- SafeSense® – autonomous, safe bed exit with individual call notification.
- 3-stop strategy – situation- and need-based care, provision and support
- Underbed clearance for mobile patient lifters at any position
- Quality German made construction



Specifications:

External Dimensions	2080 x 1060 mm
Sleeping Surface	870/900 x 2000 mm
Backrest Adjustment	70° ± 5° (retraction by 120 mm)
Upper Backrest Adjustment	30° ± 5° (retraction by 50 mm)
Head/ Foot Swivel	17° / 17°
Height	230 – 800 mm
Bed Weight	115 kg (depending on model)
SWL	270 kg