

Secura 1 Mental Health Bed

Product Code: 7-SEC1

The Secura 1 is a robust bed for demanding environments, providing a sleeping surface which can reduce harm to both the patient and carers. There are no electrics included, hence creating a safe solution for your facility, with a SWL of 200kg and durable framework. Suitable for your medium-high care patients.




Product Features:

- Australian Made and Owned
- No electrics
- Fixed height solid frame
- Secure bolt on head and footboards
- Strong robust construction



Enclosed
Mattress Supports



Anti Tamper
Screws



Natural Oak
Laminated Endboards



Ironstone
Powdercoat

Specifications:

Size	Single	King Single
Overall Dimensions	920 x 2070 mm	1090 x 2070 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm
Mattress Base Height	450 mm	450 mm
SWL	200 kg	200 kg

Optional Features – Added at time of manufacture

- Custom powder coat Colours
- Custom laminate Colours
- Built in 200mm extension
- Extra length
- Floor mounting plates



Standard Laminate Natural Oak
Standard Powdercoat Ironstone



Custom Laminate Oyster Grey
Custom Powdercoat APO Grey