

# Alrick 2400 Acute Care Bed

Product Code: 7-24

The 2400 Acute Care Bed is engineered with a robust design to ensure safety and stability. It boasts an impressive safe working load capacity of 350kg to 400kg.



## Specifications:

Size	King Single	Extra Wide	Double	Queen Double
External Dimensions	1165 x 2210 mm	1295 x 2210 mm	1465 x 2210 mm	1595 x 2210 mm
Sleeping Surface	1070 x 2000 mm	1200 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm
Hi-Lo Position	390 - 830 mm	390 - 830 mm	390 - 830 mm	390 - 830 mm
Underbed Clearance	160 mm	160 mm	160 mm	160 mm
Back / Knee Angle	77° / 40°	77° / 40°	77° / 40°	77° / 40°
TR / ATR Angle	16° / 16°	16° / 16°	16° / 16°	16° / 16°
SWL	400 kg	400 kg	350 kg	350 kg

# Ensuring safety and stability in an Acute Care environment.



## Product Features:

- Australian Made and Owned
- Vertical Lift
- Electric Auto-Regression Backrest
- CPR Backrest
- Electric Kneebreak
- Manual Vascular Support
- 350 - 400kg Safe Working Load
- Back Lit Handset
- Battery Back Up
- Under Bed Light
- Ironstone Powder Coated End Boards with White Laminate Infills
- APO Grey and Ironstone Powder Coated Frame



Forward & Reverse  
**Trendelenburg**




**Chair Position**  
Function



**150mm Central**  
Locking Castors



**Bluetooth**  
**Connectivity**  
\*\*\*handset must be  
ordered separately\*\*\*



## Optional Features - Added at time of manufacture

- Electric Vascular Support
- Pull Out Extension
- Attendant Control
- Split Platform Upgrade
- USB Charging Ports
- Custom Powder Coat Colour
- Custom End Board Colours
- Custom Length of Sleeping Surface (2m+)

## Accessories - Added after manufacture

- Mattresses
- Horizontal 400HR Rails (pair)
- Padded Cotside Protectors (pair)
- Self Help Pole
- IV Pole