

Alrick Luxa Homecare Bed

Product Code: 7-LUXA

The Luxa Homecare Bed is a premium homecare bed that is designed to give the user ultimate comfort and care, while maintaining dignity and style in the bedroom.

This bed is packed with premium features such as electric auto-regression backrest, manual vascular support, a high height of 725mm, a low bed height of 150mm and optional bluetooth handsets.



Product Features:

- Australian Made and Owned
- Floorline Bed
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- 250kg SWL
- Back Lit Handset
- Battery Back Up
- Ironstone Powder Coated Frame
- Full Surround Design



Forward & Reverse Trendelenburg



Chair Position Function






75mm Twin Wheel Braking Castors



Bluetooth Connectivity
handset must be ordered separately

Optional Features - Added at time of manufacture

- Dual Locking Castors
- Bluetooth Handset
- Under Bed Light
- USB Charging Ports
- Timber Surround
- Custom Powder Coat Colours
- Custom End Board Styles and Colours
- Custom Sleeping Surface Length (2m+)
- Partner Bed Set Up



Stacked Styles & Colours

Custom styles & colours available upon request.



Upholstered

Plain Padded Headboard & Plain Full Surround



Upholstered

Jersey Padded Headboard & Plain Full Surround




Foot section of surround available in Low or High Profile.



Natural Oak
Laminate



Ocean
Upholstered



Ash
Upholstered



Accessories - Added after manufacture

- Mattresses
- Luxa Horizontal 400HR Bed Rails (pair)
- Padded Cotside Protector (pair)
- Luxa Clamp on Support Rail
- Luxa 200mm Bed Extension
- 200mm Foam Bolster (order with Bed Extension)
- Emergency Bed Power Unit
- IV Pole
- EN T-Bar buffer
- Reading Lamp



Specifications:

Size	Single	King Single	Double	Queen Double
Overall Dimensions	1000 x 2160 mm	1170 x 2160 mm	1470 x 2160 mm	1600 x 2160 mm
Sleeping Surface	900 x 2030 mm	1070 x 2030 mm	1370 x 2030 mm	1500 x 2030 mm
Hi - Lo Position	150 - 725 mm	150 - 725 mm	170 - 725 mm	170 - 725 mm
TR / ATR Angle	62° / 35°	62° / 35°	62° / 35°	62° / 35°
Back / Knee Angle	15° / 15°	15° / 15°	15° / 15°	15° / 15°
SWL	250 kg	250 kg	250 kg	250 kg