

Alrick 2300 Tilt Care Bed

Product Code: 7-23T

The Alrick 2300 Tilt Care Bed is the most versatile Hi-Low and tilt bed available for the healthcare industry. Its robust frame supports up to 300kg SWL, accommodating users with diverse weight requirements. Available in sizes ranging from single to king double width, this bed is equally suitable for both in-home and facility care settings, combining aesthetic appeal with practicality for users and carers alike.



Product Features:

- Australian Made and Owned
- Vertical Lift
- Electric Backrest
- Electric Kneebreak
- Manual Vascular Support
- 300kg SWL
- Pull Out Extension
- Back Lit Handset
- Under Bed Light
- Curved Natural Oak End Boards
- APO Grey and Ironstone Powder Coated Frame

Optional Features - *Added at time of manufacture*

- Partner Bed Set Up
- CPR Backrest
- Electric Vascular Support
- Central Locking Castors
- Bluetooth Handset
- Attendant Control
- Battery Back Up
- USB Charging Ports
- Custom Powder Coat Colour
- Custom End Board Style and Colour
- Custom Length Sleeping Surface (1.5m+)

Accessories - *Added after manufacture*

- Mattresses
- Horizontal Rails
- Padded Cotside Protector
- Clamp On Support Rail
- StandEzy Straight Handle
- StandEzy Hook Handle
- Self Help Poles
- 200mm Foam Bolster
- Emergency Bed Power Unit
- IV Pole
- Combo Grab Handle
- T-Bar Buffer
- Reading Lamp



Forward & Reverse Trendelenburg



Chair Position Function



125mm Braking & Directional Castors



Bluetooth Connectivity
handset must be ordered separately