

Tom 2 Cot

Product Code: 76750

A state-of-the-art childrens cot that is safe and provides excellent access for nursing staff.



A state-of-the-art children's cot must be safe and be able to provide excellent access for nursing staff. The electric positioning of TOM 2 children's cot responds to these and further demands in its original concept of telescopic siderails, column construction mechanism and above standard safety parameters. The Tom 2 bed fully complies with the new standard for hospital children's beds, DIN 32623:2009-11 and its construction guarantees perfect protection against children fall. The Tom 2 Cot also complies with standards: EN 60601-1, EMC 60601-1-2, EMC 60601-1-6



Product Features:

- 3 height positions for side rails, 0cm, 40cm, 80cm
- Backrest retracts when it raises, meaning zero gap
- 80kg SWL
- The control panel is equipped with an authorization security key mechanism
- The two ergonomic handrails make the release of the side rails safe.
- Thanks to the transparent bed ends, the child patient is in constant contact with their surroundings, with parents and staff always being close by.
- Column Construction ensure perfect stability, hygiene and unrestricted access to all parts of the bed.
- Nightlight
- Accessory Holder, Box and Rails
- 100mm braking castors
- Nurse control panel
- Removable head and foot ends with clear panels
- Wall bumpers
- Quality European product



Specifications:

Product Code	76750
External Dimensions	830 x 1510 mm
Mattress Dimensions	700/100 x 1370 mm
Backrest Angle	41°
Height Range	580 - 830 mm
Siderail Height	800 mm
Backrest Auto-Regression	125 mm
TR / ATR	+ 15° / - 15°
Bed Weight	120 kg
SWL	80 kg

