

activeMEDICAL

# Custom Sling

Prescription Form



## Form Information

Full Name

Contact

Date

## General Sling Design

### Base Sling Design (loop fixing only)

General Purpose	Care	Junior	Walking / Standing
General Purpose Head Support	Care Head Support	Amputee	

### Base Sling Material

Mesh	Fabric	Parasilk
------	--------	----------

### Base Sling Type

	Loop	Clip	E Clip	Dual (Clip)	Dual (E-Clip)
Pivot					
Pivot Care					
Standing Transfer					
Super Soft					
Other					

### Base Sling Size

XS (Extra Small)	M (Medium)	XL (Extra Large)	Other *
S (Small)	L (Large)	XXL (Extra Extra Large)	

\* If you have selected 'other', please provide generic dimensions (A - D) plus any other relevant dimensions in the sketch box provided on the last page.

## Choose a Standard Modification


### MOD 001 - Change/ Add Torso Belt

Buckle	Neoprene/Buckle Over	Custom Length
Velcro	Neoprene/ Velcro	Custom Width

Comments


### MOD 002 - Add / Reduce Leg Length

mm	Add Leg Length
mm	Reduce Leg Length



### MOD 003 - Add / Reduce Leg Width


mm	Add Leg Width
mm	Reduce Leg Width



**MOD 004 -**  
Add / Reduce Head Support

mm Add Head Support


mm Reduce Head Support



**MOD 005 -**  
Add / Reduce Body Length

mm Add Body Length


mm Reduce Body Length



**MOD 006 -**  
Add / Reduce Body Width

mm Add Body Width

mm Reduce Body Width



NB - Body width adjustment will not affect generic external profile of sling.

**MOD 007 -**  
Add Modesty Straps

Modesty Straps keeps legs of sling together without crossing




**MOD 008 -**  
Care Sling Leg Strap

Assists in support of upper thighs for clients with low tone using care slings.




**MOD 009 -**  
Stand Up Hemi Band

Supports arms affected by hemiplegia in a standing transfer




**MOD 010 -**  
Stand Up Drop Strap

Assists clients with limited weight bearing in a standing transfer.



**Sling Sketch Box & Dimensions**



\* If you have selected 'other' in base sling size, please provide generic dimensions (A - D) plus any other relevant dimensions in the sketch box provided below.

**Additional Comments:**

# Sling Sketch and Dimensions

A=

B=

C=

D=

Disclaimer: While Custom Slings are covered by manufacturer's warranty, no responsibility is taken for their suitability and functionality.