

Secura 2 Bed

Prescription Form

Providing dignity and comfort for all mental health patient in healthcare environments



Form Information

Thank you for the opportunity to quote the Secura 2 for your facility. To ensure we quote exactly to your requirements, please tick each box accordingly & sign to approve. For any requests not listed, please note these in the Custom Requests section.

| | | |
|-----------|---------|------|
| Full Name | Contact | Date |
|-----------|---------|------|








Bed Requirements

| | Single | King Single |
|-----------------------------|---------------|----------------|
| Overall Dimensions | 920 x 2070 mm | 1090 x 2070 mm |
| Sleeping Surface Dimensions | 900 x 200 mm | 1070 x 2000 mm |
| Height | 390 - 830 mm | 390 - 830 mm |
| SWL | 300 kg | 300 kg |

Custom Dimensions



| | |
|-------------|--|
| Length (mm) | |
|-------------|--|

Optional Functions

| | | | |
|---|--|---|---|
| <p>Backrest</p>  | <p>Kneebreak</p>  | <p>Trendelenburg</p>  | <p>Battery Backup</p>  |
| <p>Central Locking Castors</p>  | <p>Key Power Lockout</p>  | <p>Bluetooth Handset</p>  | <p>Cable Length: Headboard to Power Point</p> <p>In the note section below:</p> <ol style="list-style-type: none"> Please specify height of power point from the floor. Please specify where the power point is in relation to the bed (feel free to attach photo of the room). |

Bed Frame Design

Headboard - Curved Style (25mm headboard & 32mm footboard)

| | |
|--|--|
| <p>Natural Oak</p>  | <p>Oyster Grey</p>  |
|--|--|

Comments

Please note below, any custom requests or any requirements that are not listed. Upon receiving these requests, we will assess if the request is possible and discuss these options further with you.

Approval

| | | | |
|--------------|----------|--|--|
| Full Name | | | |
| Organisation | Position | | |
| Mobile | Email | | |
| Signature | Date | | |